

Monochromators at Synchrotron SOLEIL

D. Dallé, M. Ribbens, G. Cauchon, D. Duran, S. Lecouster, F. Bertran, & All



Synchrotron Soleil, L'Orme des Merisiers, Saint-Aubin - BP 48
91192 GIF-sur-YVETTE CEDEX - FRANCE



This poster presents a general overview of the different monochromators in operation at Synchrotron SOLEIL. These monochromators cover a wide energy range which starts from the low energies with few eV to the high energies reaching 50keV in.
The high energies are selected by the help of common monochromator crystals like Channel-cuts or Double Crystal Monochromators (DCM). Some of them are water cooled, others are cryocooled by LN₂.
Some low energy monochromators are constituted of up to five interchangeable gratings while others high energy monochromators are equipped with 2 sets of crystals for example with one pair of silicon (111) crystals and another one pair of silicon (311).
The second crystal of the DCM should be a sagittal focusing crystal.
In addition, one beamline is equipped with a Quick EXAFs monochromator with adjustable amplitude of oscillations and a variable frequency.

Low Energy Monochromator




The low energies are selected by the help of plane grating like the usual type with Variable Line Spacing - Variable Groove Depth (VLS-VGD) along the grating lines and the "multilayer" grating. Some low energy monochromators are constituted of up to five interchangeable gratings.

The CASSIOPEE plane grating monochromator - 10 to 1500eV

The working principle of the monochromator is sketched in the figure below. One can choose the incident angle α of the white beam on the plane grating by rotating the grating around an horizontal axis. The β angle can be chosen by rotating the plane mirror, above the grating, around a horizontal axis located below the grating plane, to ensure a constant height of the outgoing beam. The CASSIOPEE monochromator contains three different gratings with 400 l/mm, 800 l/mm and 1600 l/mm respectively.

The CASSIOPEE monochromator is designed to work in the "modified Petersen mode", which means that α and β are always chosen to fulfill


$$\frac{\sin \beta}{\sin \alpha} = \epsilon \approx 0.2 \quad \text{to ensure an optimal resolution on the whole photon energy.}$$






Cassiopee - Figure of the C factor; resolving power vs Energy

Hermes - low energy monochromator

ANTARES - low energy monochromator - Tx Translation stage



The positions and the fixations of the slide bearings have a important influence on the quality of the roll of this stage as illustrated in the figure to the side.




High Energy Monochromator


SAMBA Beamline - Quick-EXAFS monochromator - energy range \rightarrow BM / 4 to 40keV

Adjustable Oscillator from the outside of the vacuum vessel

This design concerns the possibility to adjust the oscillation amplitude without opening the vessel, in other words at most several minutes during an experiment. Then the device permits to control the angular position of the crystals on the fly.
The oscillator moves a channel-cut from a frequency of 0.1 hertz on a range $\pm 2^\circ$ to a frequency of 50 hertz on a range $\pm 0.1^\circ$.




Single Quick Scan Nickel foil compared to step by step scan at various frequencies of data collection





Data quality is stable with frequency and no systematic error is introduced by the device as compared to a spectrum obtained step by step with another monochromator on the same Ni foil on the same beamline.
The only drawback is a poorer resolution beyond 10Hz inherent to the frequency cutoff of the employed detectors (ionization chambers).

Proxima 1 Beamline - Channel-Cut monochromator

- energy range \rightarrow ID / 5 to 15keV




PROXIMA 1 Monochromator results




CRISTAL Beamline - double DCM monochromator

- energy range \rightarrow ID / 4 to 30keV



Translation stage to switch from Si 1,1,1 to Si 3,1,1



Acknowledgements:

We are grateful to thank the PROXIMA, CASSIOPEE, SAMBA and CRISTAL teams for giving us information and recorded data.