

Bending Magnet and 3-Pole Wiggler Frontend Design



Prepared by: Muhammad Aftab Hussain

October, 2014



MEDSI 2014

Mechanical Engineering Design of
Synchrotron Radiation Equipment
and Instrumentation



BROOKHAVEN
NATIONAL LABORATORY


Outline

- Introduction
- Frontend Configuration & Interfaces
- Frontend Overall Design Approach and Challenges
- Frontend Components Design
- Summary and Conclusion

Introduction

Frontend serves the following purpose:

- ❑ It transports the required synchrotron radiation from the source to the first optics enclosure.
- ❑ It protects the storage ring vacuum.
- ❑ It removes the undesired heat from the synchrotron beam.
- ❑ It collimates the Bremsstrahlung radiation to a safe level.




Beamline Source Parameters

Sources for the NSLS-II BM port beamlines:

- Source Type: 3PW and BM***
- Nominal magnetic field at source location: 3PW (1.2 T) & BM (0.40 T)***
- Usable photon energy range: 1-31 keV***
- Maximum fan power (500mA current): 3PW (250 W) & BM (69 W)***
- Maximum Power Density: 3PW (365 W/mrad²) & BM (115 W/mrad²)***
- Beamline Horizontal Acceptance Fan: 0.2-3 mrad***
- Beamline Vertical Acceptance Fan: 0.15-0.6 mrad***

Frontend Configuration & Interfaces



Frontend Overall Design Approach & Challenges


Challenges

- ❑ Full passive protection to the frontend components.
- ❑ Geometric Envelope as opposed to Active Interlock Envelope for full passive protection.

Design Approach

- ❑ Ray tracing to locate and size the components.
- ❑ Multiple flange absorbers for the progressive trimming of synchrotron beam.
- ❑ Fixed mask and flange absorbers are double sided conflate flanges. (BNL Patent Pending)
- ❑ Glidcop AL-15 and Copper Chromium Zirconium as a flange material.

Frontend Configuration




MEDSI 2014

Mechanical Engineering Design of
Synchrotron Radiation Equipment
and Instrumentation

Frontend Ray Tracing

Horizontal Synchrotron Ray Tracing



Glidcop VS Cu-Cr-Zr

Material Properties

- *Thermal Conductivity (RT):*
Glidcop Al25, Al15: 344 - 365 W/(m.K)
Cu-Cr-Zr: 314 - 335 W/(m.K)
- *Elastic Modulus:*
Glidcop Al15, Al25: 130 GPa
Cu-Cr-Zr: 123 GPa
- *0.2 % Yield Strength, (RT, Cold Worked):*
Glidcop Al15, Al25: 470 - 580 MPa
Cu-Cr-Zr: 350 - 550 Mpa
- *Coefficient of Thermal Expansion:*
Glidcop Al15, Al25: 16.6 $\mu\text{m}/\text{K}$
Cu-Cr-Zr: 17.0 $\mu\text{m}/\text{K}$

- *Cu-Cr-Zr (C18150) is 1/4th the price of Glidcop AL-15.*
- *Cu-Cr-Zr is readily available in different forms and sizes from many suppliers.*
- *Cu-Cr-Zr loses its strength rapidly if exposed to sustained temperatures > 500°C*
- *Glidcop is the choice if brazing is required.*


Ref: Li M. and Zinkle S. J. (2012) Physical and Mechanical Properties of Copper and Copper alloys, Comprehensive Nuclear Materials, Vol. 4, pp 667-690

(Tests on Cu-Cr-Zr Conflate Flanges and Flare Fitting Connections)

Ref: Sushil Sharma, et al. "A Novel Design of High Power Masks & Slits", MEDSI2014


Fixed Mask / Flange Absorber Design

Fixed Mask



Lead Collimator

- Fixed mask defines the beam size and shadow collimator aperture.
- Fixed mask is water cooled and interlocked to the machine EPS.
- External fins are added to fixed mask to provide natural convection and equipment protection from loss of water flow.



Flaretite Fitting

Double Sided
CF Flange


Fins

Fixed Mask /Lead Collimator (BC1) Assembly


Fixed Mask Thermal Analysis

Analysis Parameters

- Distance from Source ~ 7.9m
- Total Power (3PW & BM) = 319 W
- Assumed heat transfer coefficient of 10 W/m².K for natural convection.
- Assumed film coefficient of 0.4 W/cm².°C for a flow rate of 1GPM in Ø0.375" cooling channel.



Temperature Contour Plot, Tmax ~ 87°C




Temperature Contour Plot, Tmax ~ 93°C


Fixed Mask Thermal Analysis (Natural Convection Only)

Analysis Parameters

- Distance from Source ~ 7.9m
- Total Power (3PW & BM) = 319 W
- Assumed heat transfer coefficient of 10 W/m².K for natural convection.
- Typically the heat transfer coefficient value for natural convection ranges from 5-25 W/m².K




Temperature Contour Plot, T_{max} ~ 298°C



Temperature Contour Plot, T_{max} ~ 305°C

Beryllium Window Design


Be Window



Glidcop AL-15
Double Sided CF Flange

Be Window/Ion Pump Assembly

- *Water cooled Be Window.*
- *Provide isolation to storage ring vacuum from beamline vacuum.*
- *Window thickness (250 μm).*
- *Window aperture (10mm x 92mm).*
- *Polished to RA < 100 nm RMS on both sides.*
- *Diffusion bonded assembly bakeable to 450°C.*
- *Capable of sustaining > 2.4 bar pressure differential.*




New Design

304L SST
CF Flange

OFE Cu
Brazed to
SST Flange

304L SST
CF Flange



Traditional Design




MEDSI 2014

Mechanical Engineering Design of
Synchrotron Radiation Equipment
and Instrumentation



BROOKHAVEN
NATIONAL LABORATORY

Be Window Test Results



Prototype Be Window

- Leak check before and after the diffusion bonding at room temperature.
- Leak check and pressure test was performed after repeated thermal cycling up to 150°C.
- Pressure test was performed to check the window strength.

NSLS-II NexGen Beamline Front-End Beryllium Window Leak and Pressure Tests			
E. Hu (5-29-2014)			
Be-Window Leak Check and Pressure Test	Valve-1 /Valve-2	Pressure-1/Pressure -2	Leak at Be-window or knife edge (side 1 or 2)
(2 identical independent sets at both sides of Be)	(open or close)	(mTorr vs 1.5 atm)	< 2E-10 (Torr-Liter/sec)
(temperature rise at rate of 50-deg C per hour)			
Leak/P-test at Room Temp using helium			
pump side-1 to LV, feed He to side-2 (1.5 atm)	open / close	mTorr / 1.5 atm	< 2x10 ⁻¹⁰
pump side-2 to LV, feed He to side-1 (1.5 atm)	close / open	1.5 atm / mTorr	< 2x10 ⁻¹⁰
Leak/P-test at 70-deg C using helium			
pump side-1 to LV, feed He to side-2 (1.5 atm)	open / close	mTorr / 1.5 atm	< 2x10 ⁻¹⁰
pump side-2 to LV, feed He to side-1 (1.5 atm)	close / open	1.5 atm / mTorr	< 2x10 ⁻¹⁰
(temperature rise at rate of 50-deg C per hour)			
Leak/P-test at 120-deg C using helium			
pump side-1 to LV, feed He to side-2 (1.5 atm)	open / close	mTorr / 1.5 atm	< 2x10 ⁻¹⁰
pump side-2 to LV, feed He to side-1 (1.5 atm)	close / open	1.5 atm / mTorr	< 2x10 ⁻¹⁰




Test Setup







X-Y Slits Design

X-Y Slits



Thermally Stable
Invar Stand

- *Water cooled , normal incidence design.*
- *Allow full independent adjustment of all four blades.*
- *Maximum Opening (50 mm H x 10 mm V)*
- *Minimum Aperture Size (-10 mm, fully closed, overlapped)*
- *Allow selection of any part of the beam.*
- *Position Resolution < 1 μm*
- *Aperture Stability (1 μm or better over any 24-hour period)*
- *Repeatability (5 μm)*
- *Beam intensity monitoring.*



X-Y Slits Assembly

Front View

Back View



MEDSI 2014


Mechanical Engineering Design of
Synchrotron Radiation Equipment
and Instrumentation




BROOKHAVEN
NATIONAL LABORATORY

Photon Shutter Design


Photon Shutter



- *Water cooled , normal incidence design.*
- *Stops synchrotron beam whenever safety shutters are closed.*
- *Integrated into the machine PPS.*



Photon Shutter
Open Position



Photon Shutter
Closed Position

Photon Shutter/Bremsstrahlung Collimator Assembly

Summary & Conclusion

- *Full passive protection.*
 - Any beam deviation is acceptable for frontend components.
 - Progressive trimming of fan using flange absorbers.
 - Flange absorbers shadow lead collimators.
 - Components protection in absence of water flow by natural convection.
- *Cost effective design.*
 - No welding or brazing required for mask, slits, absorbers & photon shutter.
 - Cu-Cr-Zr is 1/4th the price of Glidcop AL-15.
 - Cu-Cr-Zr is available in different forms (bars, sheets & plates) and sizes.
 - Standardized components for mass production.
- Glidcop AL-15 double sided CF flange is needed for Be window.
- Design work in progress for five NSLS-II BM port beamlines
 - Two frontends will be installed in 2015.

Acknowledgement

- ❖ Dick Hseuh
- ❖ Eugene Hu
- ❖ Frank DePaola
- ❖ Frank Lincoln
- ❖ John Tuozzolo
- ❖ Jeff Keister
- ❖ Paul Northrup
- ❖ Sushil Sharma

Thank you for your attention!